



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SUCNR1 Antibody [#G748]

Catalog Number: G748

Concentration: 1mg/ml

Swiss-Prot No.: Q9BXA5

Other Names: Succinate receptor 1; G-protein coupled receptor 91; P2Y purinoceptor 1-like; GPR91

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SUCNR1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SUCNR1 antibody detects endogenous levels of total SUCNR1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Wittenberger T., J. Mol. Biol. 307:799-813(2001).

Zhang W., Submitted (MAR-2000) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

