



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TMEM185A Antibody [#G749]

Catalog Number: G749

Concentration: 1mg/ml

Swiss-Prot No.: Q8NFB2

Other Names: Transmembrane protein 185A; Protein FAM11A; CXorf13, FAM11A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TMEM185A.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TMEM185A antibody detects endogenous levels of total TMEM185A protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Shaw M.A., Eur. J. Hum. Genet. 10:767-772(2002).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Ota T., Nat. Genet. 36:40-45(2004).



Western blot analysis of extracts from K562 cells, using TMEM185A antibody.



Immunofluorescence analysis of MCF-7 cells, using TMEM185A antibody.