



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## TAS2R48 Antibody [#G758]

**Catalog Number:** G758

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P59542

**Other Names:** Taste receptor type 2 member 48; T2R48

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human TAS2R48.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TAS2R48 antibody detects endogenous levels of total TAS2R48 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Bufe B., Nat. Genet. 32:397-401(2002).

Fischer A., Mol. Biol. Evol. 22:432-436(2005).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

