



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

TAS2R49 Antibody [#G759]

Catalog Number: G759

Concentration: 1mg/ml

Swiss-Prot No.: P59543

Other Names: Taste receptor type 2 member 49; T2R49; T2R56

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TAS2R49.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TAS2R49 antibody detects endogenous levels of total TAS2R49 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Bufe B., Nat. Genet. 32:397-401(2002).

Conte C., Cytogenet. Genome Res. 98:45-53(2002).

Fischer A., Mol. Biol. Evol. 22:432-436(2005).

