



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TAS2R8 Antibody [#G768]

Catalog Number: G768

Concentration: 1mg/ml

Swiss-Prot No.: Q9NYW2

Other Names: T2R8; TA2R8; taste receptor family B member 5; taste receptor type 2 member 8; taste receptor, family B, member 5; taste receptor, type 2, member 8; TRB5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

TAS2R8.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TAS2R8 antibody detects endogenous levels of total TAS2R8 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Adler E., Cell 100:693-702(2000).

Fischer A., Mol. Biol. Evol. 22:432-436(2005).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

