



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AVPR1B Antibody [#G787]

Catalog Number: G787

Concentration: 1mg/ml

Swiss-Prot No.: P47901

Other Names: antidiuretic hormone receptor 1b; arginine vasopressin receptor 1B; arginine vasopressin receptor 3; AVPR V1b; AVPR V3; AVPR3; pituitary vasopressin receptor 3; V1bR; vasopressin V1b receptor; vasopressin V3 receptor; VPR3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human AVPR1B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AVPR1B antibody detects endogenous levels of total AVPR1B protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Sugimoto T., J. Biol. Chem. 269:27088-27092(1994).
de Keyzer Y., FEBS Lett. 356:215-220(1994).
North W.G., Peptides 18:985-993(1997).

