



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

VN1R2 Antibody [#G792]

Catalog Number: G792

Concentration: 1mg/ml

Swiss-Prot No.: Q8NFZ6

Other Names: Vomeronasal type-1 receptor 2; V1r-like receptor 2; hGPCR25; V1RL2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human VN1R2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: VN1R2 antibody detects endogenous levels of total VN1R2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Rodriguez I., Curr. Biol. 12:R409-R411(2002).

Zhang J., Proc. Natl. Acad. Sci. U.S.A. 100:8337-8341(2003).

Takeda S., FEBS Lett. 520:97-101(2002).

