



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

VN1R4 Antibody [#G794]

Catalog Number: G794

Concentration: 1mg/ml

Swiss-Prot No.: Q7Z5H5

Other Names: Vomeronasal type-1 receptor 4; V1r-like receptor 4; hGPCR27; V1RL4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human VN1R4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: VN1R4 antibody detects endogenous levels of total VN1R4 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Rodriguez I., *Curr. Biol.* 12:R409-R411(2002).

Zhang J., *Proc. Natl. Acad. Sci. U.S.A.* 100:8337-8341(2003).

Takeda S., *FEBS Lett.* 520:97-101(2002).



Western blot analysis of extracts from COS-7 cells, using VN1R4 antibody.



Immunofluorescence analysis of MCF-7 cells, using VN1R4 antibody.