



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

OR10A5 Antibody [#G819]

Catalog Number: G819

Concentration: 1mg/ml

Swiss-Prot No.: Q9H207

Other Names: Olfactory receptor 10A5; Olfactory receptor-like protein JCG6; Olfactory receptor 11-403; OR11-403; HP3; R10A1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR10A5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR10A5 antibody detects endogenous levels of total OR10A5 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:20000

References:

Gaudin J.-C., Chem. Senses 26:1157-1166(2001).

Lane R.P., Proc. Natl. Acad. Sci. U.S.A. 98:7390-7395(2001).

Buettner J.A., Genomics 53:56-68(1998).

