



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

OR10G7 Antibody [#G825]

Catalog Number: G825

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGN6

Other Names: O10G7; olfactory receptor 10G7;
OR11-283

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR10G7.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: OR10G7 antibody detects endogenous levels of total OR10G7 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.



Western blot analysis of extracts from HepG2 cells, using OR10G7 antibody.



Immunofluorescence analysis of A549 cells, using OR10G7 antibody.