



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

OR10H1 Antibody [#G827]

Catalog Number: G827

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y4A9

Other Names: O10H1; olfactory receptor 10H1; olfactory receptor, family 10, subfamily H, member 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OR10H1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR10H1 antibody detects endogenous levels of total OR10H1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Grimwood J., Nature 428:529-535(2004).

Fuchs T., Genomics 80:295-302(2002).

