



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

OR10J3 Antibody [#G828]

Catalog Number: G828

Concentration: 1mg/ml

Swiss-Prot No.: Q5JRS4

Other Names: Olfactory receptor 10J3; OR10J3;
OR10J3P

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OR10J3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR10J3 antibody detects endogenous levels of total OR10J3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Gregory S.G., Nature 441:315-321(2006).

