



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## OR10R2 Antibody [#G830]

**Catalog Number:** G830

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q8NGX6

**Other Names:** O10R2; olfactory receptor 10R2; olfactory receptor, family 10, subfamily R, member 2; OR1-8; OR10R2Q

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human OR10R2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR10R2 antibody detects endogenous levels of total OR10R2 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Gregory S.G., Nature 441:315-321(2006).

Fuchs T., Genomics 80:295-302(2002).

