



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

OR51H1 Antibody [#G848]

Catalog Number: G848

Concentration: 1mg/ml

Swiss-Prot No.: Q8NH63

Other Names: Putative olfactory receptor 51H1; Olfactory receptor OR11-25; OR51H1P; OR51H1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR51H1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR51H1 antibody detects endogenous levels of total OR51H1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:7205-7205(2004).

