



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## OR52A5 Antibody [#G851]

**Catalog Number:** G851

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9H2C5

**Other Names:** Olfactory receptor 52A5; Olfactory receptor OR11-33; Odorant receptor HOR3'beta5

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human OR52A5.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR52A5 antibody detects endogenous levels of total OR52A5 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Bulger M., Proc. Natl. Acad. Sci. U.S.A. 97:14560-14565(2000).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

