



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## OR2F2 Antibody [#G893]

**Catalog Number:** G893

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O95006

**Other Names:** olfactory receptor, family 2, subfamily f, member 2; or2f2; or7-1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human OR2F2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR2F2 antibody detects endogenous levels of total OR2F2 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Hillier L.W., Nature 424:157-164(2003).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

