



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ADAM 17 (Cleaved-Arg215) Antibody [#L0130]

**Catalog Number:** L0130

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P78536

**Other Names:** a disintegrin and metalloproteinase domain 17; ADA17; ADAM 17; CD156b antigen; CSVP; EC 3.4.24.86; snake venom-like protease; TACE; TNF-alpha convertase; TNF-alpha-converting enzyme

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from human ADAM 17.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ADAM 17 (Cleaved-Arg215) antibody detects endogenous levels of fragment of activated ADAM 17 resulting from cleavage adjacent to Arg215.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Moss M.L., Nature 385:733-736(1997).

Black R.A., Nature 385:729-733(1997).

Patel I.R., J. Immunol. 160:4570-4579(1998).

