



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CASP3 (p17, Cleaved-Asp175) Antibody [#L0153]

Catalog Number: L0153

Concentration: 1mg/ml

Swiss-Prot No.: P42574

Other Names: apopain; CASP-3; Caspase-3; CPP-32; CPP32; cysteine protease CPP32; EC 3.4.22.56; ICE3; SCA-1; SREBP cleavage activity 1; yama

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human Caspase 3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: CASP3 (p17, Cleaved-Asp175) antibody detects endogenous levels of fragment of activated Caspase 3 resulting from cleavage adjacent to Asp175.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Fernandes-Alnemri T., J. Biol. Chem. 269:30761-30764(1994).

Tewari M., Cell 81:801-809(1995).

Pelletier M., Biochem. Biophys. Res. Commun. 316:93-99(2004).



For Research Use Only