



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PRDM12 Antibody [#C11694]

Catalog Number: C11694

Concentration: 1mg/ml

Swiss-Prot No.: Q9H4Q4

Other Names: PR domain-containing protein 12, PR domain zinc finger protein 12, PRDM12, PFM9

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

PRDM12.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PRDM12 antibody detects endogenous levels of total PRDM12 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Yang X.-H., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

