



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## FKRP Antibody [#C11711]

**Catalog Number:** C11711

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9H9S5

**Other Names:** Fukutin-related protein, EC 2.-.-.-, FKRP

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human FKRP.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** FKRP antibody detects endogenous levels of total FKRP protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Brockington M., Am. J. Hum. Genet. 69:1198-1209(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Brockington M., Hum. Mol. Genet. 10:2851-2859(2001).

