



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LASS4 Antibody [#C11714]

Catalog Number: C11714

Concentration: 1mg/ml

Swiss-Prot No.: Q9HA82

Other Names: LAG1 longevity assurance homolog 4, LASS4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human LASS4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LASS4 antibody detects endogenous levels of total LASS4 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:1000

References:


Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using LASS4 antibody.



Western blot analysis of extracts from COLO cells, Jurkat cells and HUVEC cells, using LASS4 antibody.