



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CLK4 Antibody [#C11717]

Catalog Number: C11717

Concentration: 1mg/ml

Swiss-Prot No.: Q9HAZ1

Other Names: Dual specificity protein kinase CLK4, EC 2.7.12.1, CDC-like kinase 4, CLK4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CLK4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CLK4 antibody detects endogenous levels of total CLK4 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Schultz J., Genomics 71:368-370(2001).

McNally T., Protein Sci. 12:1562-1566(2003).

Kantham L., Biochem. Biophys. Res. Commun. 304:125-129(2003).

