



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARNT2 Antibody [#C11723]

Catalog Number: C11723

Concentration: 1mg/ml

Swiss-Prot No.: Q9HBZ2

Other Names: ARNT protein 2, Aryl hydrocarbon receptor nuclear translocator 2, ARNT2, KIAA0307

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ARNT2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARNT2 antibody detects endogenous levels of total ARNT2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Barrow L.L., Teratology 66:85-90(2002).

Nagase T., DNA Res. 4:141-150(1997).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

