



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CKI-γ1 Antibody [#C11728]

Catalog Number: C11728

Concentration: 1mg/ml

Swiss-Prot No.: Q9HCP0

Other Names: EC 2.7.11.1, CKI-gamma 1, Casein kinase I isoform gamma-1, CSNK1G1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human CKI-γ1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CKI-γ1 antibody detects endogenous levels of total CKI-γ1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kusuda J., Cytogenet. Cell Genet. 90:298-302(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Hillman R.T., Genome Biol. 5:RESEARCH008.1-RESEARCH008.16(2004).

