



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CREBZF Antibody [#C11755]

Catalog Number: C11755

Concentration: 1mg/ml

Swiss-Prot No.: Q9NS37

Other Names: CREB/ATF bZIP transcription factor, Host cell factor-binding transcription factor Zhangfei, HCF-binding transcription factor Zhangfei, CREBZF, ZF

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CREBZF.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CREBZF antibody detects endogenous levels of total CREBZF protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Lu R., Nucleic Acids Res. 28:2446-2454(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Misra V., J. Biol. Chem. 280:15257-15266(2005).



Western blot analysis of extracts from RAW264.7 cells, using CREBZF antibody.



Immunofluorescence analysis of HepG2 cells, using CREBZF antibody.