



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARTS-1 Antibody [#C11786]

Catalog Number: C11786

Concentration: 1mg/ml

Swiss-Prot No.: Q9NZ08

Other Names: A- LAP, APPILS, ARTS-1, Adipocyte-derived leucine aminopeptidase precursor, Aminopeptidase PILS, EC 3.4.11.-, PILS-AP, Puromycin-insensitive leucyl- specific aminopeptidase, Type 1 tumor necrosis factor receptor shedding aminopeptidase regulator

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human ARTS-1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARTS-1 antibody detects endogenous levels of total ARTS-1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Hattori A., J. Biochem. 130:235-241(2001).

Nakajima D., DNA Res. 9:99-106(2002).

Clark H.F., Genome Res. 13:2265-2270(2003).

