



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## HMG20B Antibody [#C11798]

**Catalog Number:** C11798

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9P0W2

**Other Names:** HMGX1, high-mobility group 20A

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human HMG20B.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** HMG20B antibody detects endogenous levels of total HMG20B protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Sumoy L., Cytogenet. Cell Genet. 88:62-67(2000).

Marmorstein L.Y., Cell 104:247-257(2001).

Lee Y.M., Biochim. Biophys. Acta 1574:269-276(2002).

