



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Collagen XI α 1 Antibody [#C12214]

Catalog Number: C12214

Concentration: 1mg/ml

Swiss-Prot No.: P12107

Other Names: Collagen alpha-1(XI) chain [Precursor], COBA1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Collagen XI α 1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Collagen XI α 1 antibody detects endogenous levels of total Collagen XI α 1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Yoshioka H., J. Biol. Chem. 265:6423-6426(1990).

Annunen S., Am. J. Hum. Genet. 65:974-983(1999).

Zhidkova N.I., J. Biol. Chem. 270:9486-9493(1995).

