



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cofilin Antibody [#C12219]

Catalog Number: C12219

Concentration: 1mg/ml

Swiss-Prot No.: P23528/Q9Y281

Other Names: CFL2, COF2, Cofilin 2, Cofilin, muscle isoform

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cofilin.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cofilin antibody detects endogenous levels of total Cofilin protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Thirion C., Eur. J. Biochem. 268:3473-3482(2001).

Zhang Y., Mol. Cell. Proteomics 4:1240-1250(2005).

Olsen J.V., Cell 127:635-648(2006).

