



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Collagen XXIII α 1 Antibody [#C12226]

Catalog Number: C12226

Concentration: 1mg/ml

Swiss-Prot No.: Q86Y22

Other Names: Collagen alpha-1(XXIII) chain, COL23A1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Collagen XXIII α 1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Collagen XXIII α 1 antibody detects endogenous levels of total Collagen XXIII α 1 protein.

Reactivity: Human, Mouse, Rat

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Banyard J., J. Biol. Chem. 278:20989-20994(2003).

The MGC Project Team, Genome Res.

14:2121-2127(2004).

Bechtel S., BMC Genomics 8:399-399(2007).

