



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

COX5A Antibody [#C12237]

Catalog Number: C12237

Concentration: 1mg/ml

Swiss-Prot No.: P20674

Other Names: COX; COX-VA; COX5A; Cytochrome c oxidase polypeptide Va; cytochrome c oxidase polypeptide, mitochondrial; Cytochrome c oxidase subunit 5A, mitochondrial; cytochrome c oxidase subunit Va; mitochondrial cytochrome c oxidase subunit Va; VA

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human COX5A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: COX5A antibody detects endogenous levels of total COX5Aprotein.

Reactivity: Human; Mouse; Rat


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Rizzuto R., Gene 69:245-256(1988).

Uddin M., BMC Evol. Biol. 8:8-8(2008).

Ebert L., Submitted (MAY-2004) to the EMBL/GenBank/DDBJ databases.



For Research Use Only