



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

COX6C Antibody [#C12239]

Catalog Number: C12239

Concentration: 1mg/ml

Swiss-Prot No.: P09669

Other Names: EC 1.9.3.1, cytochrome c oxidase polypeptide VIc

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human COX6C.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: COX6C antibody detects endogenous levels of total COX6C protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Otsuka M., Nucleic Acids Res. 16:10916-10916(1988).

Ohta S., Submitted (FEB-1996) to the EMBL/GenBank/DDBJ databases.

Hofmann S., Cytogenet. Cell Genet. 83:226-227(1998).



Western blot analysis of extracts from HUVEC cells, using COX6C antibody.



Immunofluorescence analysis of HepG2 cells, using COX6C antibody.