



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Collagen XIX α 1 Antibody [#C12224]

Catalog Number: C12224

Concentration: 1mg/ml

Swiss-Prot No.: Q14993

Other Names: Collagen alpha-1(XIX) chain, Collagen alpha-1(Y) chain, COJA1

Storage/Stability: Store at -20°C /1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Collagen XIX α 1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Collagen XIX α 1 antibody detects endogenous levels of total Collagen XIX α 1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:1000

References:

Inoguchi K., J. Biochem. 117:137-146(1995).


Khaleduzzaman M., Genomics 45:304-312(1997).

Mungall A.J., Nature 425:805-811(2003).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human thymus gland tissue using Collagen XIX α 1 antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using Collagen XIX α 1 antibody.