



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cytochrome P450 26A1 Antibody [#C12254]

Catalog Number: C12254

Concentration: 1mg/ml

Swiss-Prot No.: O43174

Other Names: Cytochrome P450 26A1, EC 1.14.-.-, Retinoic acid-metabolizing cytochrome, P450 retinoic acid-inactivating 1, P450RAI, hP450RAI, Retinoic acid 4-hydroxylase, CYP26, P450RAI1, CP26A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 26A1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cytochrome P450 26A1 antibody detects endogenous levels of total Cytochrome P450 26A1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

White J.A., J. Biol. Chem. 272:18538-18541(1997).

Sonneveld E., Cell Growth Differ. 9:629-637(1998).

ova-Griffin M.E., Biochem. Biophys. Res. Commun. 252:487-491(1998).

