



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Cytochrome P450 2A6 Antibody [#C12258]

**Catalog Number:** C12258

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P11509

**Other Names:** Cytochrome P450 2A6, EC 1.14.14.1, CYP11A6, Coumarin 7-hydroxylase, P450 IIA3, CYP2A3, P450(I), CP2A6, CYP2A3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human

Cytochrome P450 2A6.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Cytochrome P450 2A6 antibody detects endogenous levels of total Cytochrome P450 2A6 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Miles J.S., *Acids Res.* 17:2907-2917(1989).

Yamano S., *Biochemistry* 29:1322-1329(1990).

Saito S., *J. Hum. Genet.* 48:249-270(2003).

