



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cytochrome P450 2A13 Antibody [#C12260]

Catalog Number: C12260

Concentration: 1mg/ml

Swiss-Prot No.: Q16696

Other Names: Cytochrome P450 2A13, EC 1.14.14.1, CYP11A13, CP2AD

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 2A13.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cytochrome P450 2A13 antibody detects endogenous levels of total Cytochrome P450 2A13 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Fernandez-Salguero P., Am. J. Hum. Genet. 57:651-660(1995).

Su T., cer Res. 60:5074-5079(2000).

Cauffiez C., m. Biophys. Res. Commun. 317:662-669(2004).

