



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Cytochrome P450 2B6 Antibody [#C12261]

**Catalog Number:** C12261

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P20813

**Other Names:** CP2B6, CPB1, CYP2B-1, CYP2B6, CYPIIB1, CYPIIB6, Cytochrome P450 2B1, Cytochrome P450 2B6, EC 1.14.14.1, P450 IIB1, P450-B, P450-LM2, P450-PB1, P450-PB2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

Cytochrome P450 2B6.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Cytochrome P450 2B6 antibody detects endogenous levels of total Cytochrome P450 2B6 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Yamano S., Biochemistry 28:7340-7348(1989).

Grimwood J., Nature 428:529-535(2004).

Thum T., Lancet 355:979-983(2000).

