



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cytochrome P450 2C8/9/18/19 Antibody [#C12264]

Catalog Number: C12264

Concentration: 1mg/ml

Swiss-Prot No.: P33260/P10632/P11712/P33261

Other Names: CP2CI, CPCI, CYPIIC18, EC 1.14.14.1, P450-6B/29C, cytochrome P450 2C18

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 2C8/9/18/19.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cytochrome P450 2C8/9/18/19 antibody detects endogenous levels of total Cytochrome P450 2C8/9/18/19 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Romkes M., *Biochemistry* 30:3247-3255(1991).

de Morais S.M., *Biochem. Biophys. Res. Commun.* 194:194-201(1993).

Deloukas P., *Nature* 429:375-381(2004).

