



Web: www.anbio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cytochrome P450 2W1 Antibody [#C12273]

Catalog Number: C12273

Concentration: 1mg/ml

Swiss-Prot No.: Q8TAV3

Other Names: Cytochrome P450 2W1, EC 1.14.14.-, CYPIIW1, CP2W1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 2W1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: Cytochrome P450 2W1 antibody detects endogenous levels of total Cytochrome P450 2W1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:40000


References:

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using Cytochrome P450 2W1 antibody.



Peptide - +

Immunofluorescence analysis of HuvEc cells, using Cytochrome P450 2W1 antibody.