



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cytochrome P450 39A1 Antibody [#C12275]

Catalog Number: C12275

Concentration: 1mg/ml

Swiss-Prot No.: Q9NYL5

Other Names: Cytochrome P450 39A1, EC 1.14.13.99, 24-hydroxycholesterol 7-alpha-hydroxylase, Oxysterol 7-alpha-hydroxylase, Hcyp39a1, CP39A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 39A1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cytochrome P450 39A1 antibody detects endogenous levels of total Cytochrome P450 39A1 protein.


Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000


References:

Li-Hawkins J., . Biol. Chem. 275:16543-16549(2000).

The MGC Project Team, ome Res. 14:2121-2127(2004).



Western blot analysis of extracts from LOVO cells and HUVEC cells, using Cytochrome P450 39A1 antibody.



Immunofluorescence analysis of A549 cells, using Cytochrome P450 39A1 antibody.