



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cytochrome P450 4X1 Antibody [#C12290]

Catalog Number: C12290

Concentration: 1mg/ml

Swiss-Prot No.: Q8N118

Other Names: CP4X1; CYP1VX1; Cytochrome P450 4X1;
EC 1.14.14.1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 4X1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cytochrome P450 4X1 antibody detects endogenous levels of total Cytochrome P450 4X1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:40000

References:

Clark H.F., Genome Res. 13:2265-2270(2003).


Ota T., Nat. Genet. 36:40-45(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human liver carcinoma tissue using Cytochrome P450 4X1 antibody.



Western blot analysis of extracts from COS7 cells, using Cytochrome P450 4X1 antibody.