



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DNAJB4 Antibody [#C12333]

Catalog Number: C12333

Concentration: 1mg/ml

Swiss-Prot No.: Q9UDY4

Other Names: DJB4; DnaJ homolog subfamily B member 4; DNAJW; DNJB4; Heat shock 40 kDa protein 1 homolog; Heat shock protein 40 homolog; HLJ1; HSP40 homolog

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human DNAJB4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DNAJB4 antibody detects endogenous levels of total DNAJB4 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Hoe K.L., Biochim. Biophys. Acta 1383:4-8(1998).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Rush J., Nat. Biotechnol. 23:94-101(2005).

