



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FGF23 Antibody [#C12369]

Catalog Number: C12369

Concentration: 1mg/ml

Swiss-Prot No.: Q9GZV9

Other Names: Fibroblast growth factor 23, FGF-23,
Tumor-derived hypophosphatemia-inducing factor, HYPF,

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FGF23.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: FGF23 antibody detects endogenous levels of total FGF23 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 85%, Positives = 85%);
Rat (Identities = 85%, Positives = 85%)

Applications: IF: 1:100~1:500 ELISA: 1:20000

References:

Yamashita T., Biochem. Biophys. Res. Commun. 277:494-498(2000).

White K.E., Nat. Genet. 26:345-348(2000).

Shimada T., Proc. Natl. Acad. Sci. U.S.A. 98:6500-6505(2001).

