



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FRS3 Antibody [#C12384]

Catalog Number: C12384

Concentration: 1mg/ml

Swiss-Prot No.: O43559

Other Names: fibroblast growth factor receptor substrate 3; FRS2B; FRS2beta; MGC17167; SNT-2; SNT2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FRS3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FRS3 antibody detects endogenous levels of total FRS3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Xu H., J. Biol. Chem. 273:17987-17990(1998).

Mungall A.J., Nature 425:805-811(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

