



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GBP1 Antibody [#C12390]

Catalog Number: C12390

Concentration: 1mg/ml

Swiss-Prot No.: P32455

Other Names: GBP-1; Guanine nucleotide-binding protein 1; Interferon-gamma inducible protein MAG-1; Interferon-induced guanylate-binding protein 1; MAG-1; MPA1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GBP1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GBP1 antibody detects endogenous levels of total GBP1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:40000

References:


Cheng Y.-S.E., Mol. Cell. Biol. 11:4717-4725(1991).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Nantais D.E., J. Leukoc. Biol. 60:423-431(1996).



Western blot analysis of extracts from HUVEC cells, using GBP1 antibody.



Immunofluorescence analysis of MCF-7 cells, using GBP1 antibody.

For Research Use Only