



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GRID2 Antibody [#C12395]

Catalog Number: C12395

Concentration: 1mg/ml

Swiss-Prot No.: O43424

Other Names: GluR delta-2; GLURD2; glutamate receptor delta-2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GRID2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GRID2 antibody detects endogenous levels of total GRID2 protein.


Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000


References:

Hu W., Genomics 47:143-145(1998).

Sjoebloom T., Science 314:268-274(2006).



Western blot analysis of extracts from NIH/3T3 cells, using GRID2 antibody.



Immunofluorescence analysis of HuvEc cells, using GRID2 antibody.

For Research Use Only