



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IL20RB Antibody [#C12439]

Catalog Number: C12439

Concentration: 1mg/ml

Swiss-Prot No.: Q6UXL0

Other Names: DIRS1; FNDC6; I20RB; IL-20R-beta;
IL-20R2; interleukin 20 receptor beta

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IL20RB.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IL20RB antibody detects endogenous levels of total IL20RB protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Clark H.F., Genome Res. 13:2265-2270(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Zhang Z., Protein Sci. 13:2819-2824(2004).

