



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IFM3 Antibody [#C12462]

Catalog Number: C12462

Concentration: 1mg/ml

Swiss-Prot No.: Q01628

Other Names: IFM3; INIB; interferon-induced transmembrane protein 3; interferon-inducible protein 1-8U

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IFM3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IFM3 antibody detects endogenous levels of total IFM3 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:40000


References:

Lewin A.R., Eur. J. Biochem. 199:417-423(1991).

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from HepG2 cells and HeLa cells, using IFM3 antibody.



Immunofluorescence analysis of MCF-7 cells, using IFM3 antibody.