



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAPK9 Antibody [#C12684]

Catalog Number: C12684

Concentration: 1mg/ml

Swiss-Prot No.: P45984

Other Names: Mitogen-activated protein kinase 9, EC 2.7.11.24, Stress-activated protein kinase JNK2, c-Jun N-terminal kinase 2, JNK-55, MAPK9, JNK2, PRKM9

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MAPK9.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAPK9 antibody detects endogenous levels of total MAPK9 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Sluss H.K., Mol. Cell. Biol. 14:8376-8384(1994).

Kallunki T., Genes Dev. 8:2996-3007(1994).

Gupta S., EMBO J. 15:2760-2770(1996).



For Research Use Only