



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Cyclosome 1 Antibody [#C12817]

Catalog Number: C12817

Concentration: 1mg/ml

Swiss-Prot No.: Q9H1A4

Other Names: ANAPC1, ANC1, Anaphase promoting complex subunit 1, Cyclosome subunit 1, Mitotic checkpoint regulator, Protein Tsg24, TSG24

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Cyclosome 1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Cyclosome 1 antibody detects endogenous levels of total Cyclosome 1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Joergensen P.M., Gene 262:51-59(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Kraft C., EMBO J. 22:6598-6609(2003).

